




# Parity determination via Polarization

Electromagnetic waves are polarised:



The polarisation has an effect on the Compton Scattering!

# Some definitions



This and most of the following slides are based on  
B. Schlitt, RDH, et al, NIMA337 (1994) 416-426

# Klein-Nishina Formula

The scattered gamma-ray has energy


$$E'_\gamma = \frac{E_\gamma}{1 + \frac{E_\gamma}{m_0 c^2} (1 - \cos \vartheta)}$$

The differential cross section then becomes:





$$\frac{d\sigma}{d\Omega} = \frac{e^4}{2m_0^2 c^4} \frac{E'^2}{E_\gamma} \left( \frac{E'_\gamma}{E_\gamma} + \frac{E_\gamma}{E'_\gamma} - 2 \sin^2 \vartheta \cos^2 \varphi \right)$$

$r_0 = \frac{e^2}{m_0 c^2}$  is the classical electron radius. Note that  $E'_\gamma$  depends on  $\vartheta$  !

# Klein-Nishina Formula




# Compton Polarimeter



B. Schlitt, RDH, et al, NIMA337 (1994) 416-426

# Analysing Power


B. Schlitt, RDH, et al, NIMA337 (1994) 416-426



# Analysing Power II


From Klein-Nishina and the observed asymmetry

$$\varepsilon = \frac{aN_{\perp} - N_{\parallel}}{aN_{\perp} + N_{\parallel}}$$



$$\Sigma_{CE} = \frac{\frac{d\sigma}{d\Omega}(\vartheta, \varphi = 90^\circ) - \frac{d\sigma}{d\Omega}(\vartheta, \varphi = 0^\circ)}{\frac{d\sigma}{d\Omega}(\vartheta, \varphi = 90^\circ) + \frac{d\sigma}{d\Omega}(\vartheta, \varphi = 0^\circ)},$$


# Apparative asymmetry



- $^{56}\text{Co}$  source gives unpolarised gamma rays, thus we can use it to extract the apparatusive asymmetry  $a$


B. Schlitt, RDH, et al, NIMA337 (1994) 416-426

# Sample Spectra




- Singles
- Sum
- Difference at  $90^\circ$
- Difference at  $180^\circ$


# Application to $^{162}\text{Dy}$



- The concept is easily extended to different geometries, e.g. Cluster detectors.



# Gamma Ray Tracking






# AGATA

(Advanced GAMMA Tracking Array)

4 $\pi$   $\gamma$ -array for Nuclear Physics Experiments at European accelerators providing radioactive and high-intensity stable beams



## Main features of AGATA

**Efficiency:** 43% ( $M_\gamma = 1$ ) 28% ( $M_\gamma = 30$ )  
today's arrays ~10% (gain ~4) 5% (gain ~1000)

**Peak/Total:** 58% ( $M_\gamma = 1$ ) 49% ( $M_\gamma = 30$ )  
today ~55% 40%

**Angular Resolution:**  $\sim 1^\circ \rightarrow$   
FWHM (1 MeV,  $v/c=50\%$ )  $\sim 6$  keV !!!  
today  $\sim 40$  keV


**Rates:** 3 MHz ( $M_\gamma = 1$ ) 300 kHz ( $M_\gamma = 30$ )  
today 1 MHz 20 kHz



- 180 large volume 36-fold segmented Ge crystals in 60 triple-clusters
- Digital electronics and sophisticated Pulse Shape Analysis algorithms allow
- Operation of Ge detectors in position sensitive mode  $\rightarrow$   $\gamma$ -ray tracking



# Doppler Broadening




*Doppler shift*

$$E_\gamma = E_\gamma^0 \frac{\sqrt{1 - \frac{V^2}{c^2}}}{1 - \frac{V}{c} \cos \theta}$$

Broadening of detected gamma energy from

- Spread in speed of recoil  $\Delta V$
- Distribution in the direction of recoil  $\Delta \theta_N$
- Detector opening angle  $\Delta \theta_D$





Taken from D. Radford

# The Concept






Without  
Compton  
suppression  
shields



Compton  
continuum.  
=> Large  
peak to  
total ratio





With BGO  
shielding



Less solid  
angle  
coverage  
=> Big drop in  
efficiency


With highly  
segmented  
detectors

$$E'_\gamma = \frac{E_\gamma}{1 + (E_\gamma / mc^2)(1 - \cos \theta)}$$




Path of  $\gamma$ -ray reconstructed to form  
full energy event  
=> Compton continuum reduced  
=> Excellent efficiency ~50% @1MeV  
=> Greatly improved angular  
resolution (~1°) to reduce Doppler  
effects

# Ingredients of $\gamma$ -Tracking



1

Highly segmented  
HPGe detectors



2

Digital electronics  
to record and  
process segment  
signals



3


Pulse Shape Analysis  
to decompose  
recorded waves

4



Identified  
interaction

$$(x, y, z, E, t)_i$$


Reconstruction of tracks  
e.g. by evaluation of  
permutations  
of interaction points



reconstructed  $\gamma$ -rays



# AGATA





Ge crystals size:  
length 90 mm  
diameter 80 mm

**180** hexagonal crystals **3 shapes**  
60 triple-clusters all equal  
Inner radius (Ge) 22 cm  
**Amount of germanium 310 kg**  
Solid angle coverage 80 %  
Singles rate ~50 kHz  
6480 segments  
Efficiency: 40% ( $M_{\gamma}=1$ ) 25% ( $M_{\gamma}=30$ )  
Peak/Total: 65% ( $M_{\gamma}=1$ ) 50% ( $M_{\gamma}=30$ )


# Pulse Shape Analysis

- Shape of the generated signal as a function of time is dependent on the position of the interaction within the depth of the detector, which charge carrier provides the pulse, and its mobility within the detector
- Therefore, information can be gained on interaction position




# Parametric PSA


- Charge collection time for the “real charge” pulse varies according to the distance the charge carriers have to travel to the electrodes
- Rise time parameters, e.g. t<sub>30</sub> and t<sub>90</sub> can be used ( e.g. t<sub>30</sub>=T<sub>30</sub>-T<sub>10</sub> and t<sub>90</sub>=T<sub>90</sub>-T<sub>10</sub>)
- For planar detectors, with uniform electric field, gives information on the depth. For non-planar, can also give additional information, e.g. BEGe detectors give some radial information




# Parametric PSA




# Parametric PSA








Charge (arb)




Up



Left



Seg F1





Right

Time (ns)


## Parametric PSA

- Image charges are induced in the electrodes which neighbour the “sensing” electrode, due to the movement of the charge carriers
- The net image charge signal is zero because the sum of the image charges is equal and opposite in opposite electrodes
- The magnitude of the image charges in the strips (segments) either side of the sensing electrode depends on the position of interaction within the sensing electrode
- If the interaction occurs closer to one neighbouring strip than the other, the relative image charge size is larger
- In PSA, image charge size is calculated and compared, for (x-y) info



# Parametric PSA




# AGATA



# Parametric PSA



# AGATA Database




From M. Dimmock, PhD Liverpool, 2008

# Database PSA

- Generation of a database containing signals for variation interaction positions within detector
- Signals generated through experimental measurements or models
- Experimental generation is lengthy and results are unique per detector
- Vary according to detector geometry, crystallographic axis, impurity concentration
- Algorithms used to match observed signals in-situ with those from database

# Gamma-ray Tracking



# 220 MeV $^{56}\text{Fe} \rightarrow ^{197}\text{Au}$

## ATC1 + DANTE

